

# Eye Morphology TCPLA\_EYE\_002

## Purpose

To detect abnormalities in eye morphology.

## Experimental Design

- **Minimum number of animals** : 7M + 7F
- **Age at test**: Week 58
- **Sex**: We do not expect the results of this test to show sexual dimorphism

## Procedure

1. Examine the anterior of both eyes (e.g. with slit lamp) and record any abnormalities
2. Test the iris/pupil light response
3. Image abnormal eyes as a minimum or all eyes if capacity permits
4. Dilate both eyes
5. Examine the anterior and posterior of both dilated eyes (e.g. with slit lamp and ophthalmoscope) and record any abnormalities
6. Image abnormal eyes as a minimum or all eyes if capacity permits

OCT:

1. Turn on the OCT and start the database
2. Anaesthetize mouse
3. Prepare mouse eyes with drops and place contact lens (focal length 10 mm) on the right eye
4. Enter mouse data in the "Create new patient file" area and switch to the "Acquisition" window
5. Move the OCT camera to the right position and activate measurement mode
6. Place mouse collaterally to the OCT camera on the right side of a platform that is fixed in front of the OCT lens
7. Search the contact lens in the live picture of the fundus image field and place the pupil of the mouse eye in the centre of the window
8. Move the OCT camera such that OCT lens and contact lens touch each other
9. Focus the fundus picture by slightly moving up/down or forward/backward
10. Save fundus images
11. Set the "Ref.Arm" ruler such that the section of the retina is placed in the centre of the blue rectangle
12. Set the mode of measurement on "vertical, horizontal line"
13. Move the blue horizontal line in the fundus image field to the optic nerve level
14. Save images of retinal sections
15. Move the OCT camera to the left position

16. Repeat measurement procedure for the left eye

#### Scheimpflug Imaging:

1. Turn on the Pentacam and start the patient data management
2. Apply one drop 0.5% Atropine to each mouse eye for pupil dilation
3. Enter mouse data in the "Patient" group box and switch to the Scan menu
4. Activate the "1 Picture" modus in the "Image Options" area
5. Move Pentacam to the right position
6. Hold the mouse on a platform such that the vertical LED 475 nm light slit is orientated in the center of the right eye ball
7. Guarantee optimal focus by using the fine adjustment software tool in the adjustment window
8. Start imaging manually by pressing the "Start Scan" button
9. Scheimpflug images are saved automatically
10. Move Pentacam to the left position
11. Repeat measurement procedure for the left eye

## Notes

- As a minimum, all abnormalities should be imaged.
  - Where capacity permits, all mice can be imaged
- Majority of parameters can be analysed using the standard approach for assessing categorical data. To increase power for analysis purposes, where an abnormality is detected in the left, right or both eyes, the data may be combined to generate one "abnormal" category.
- Data for both eyes is recorded under one parameter to distinguish phenotypes of incomplete penetrance in individuals and if an observation for one or both eyes cannot be made, this is recorded as 'no data'. The IMPC analysis pipeline does not take into account whether an abnormality is fully penetrant or not and the same weight is given for an abnormal observations in one or both eyes. In cases where it is not possible to confirm if an abnormality is present or not, the data is not included in the statistical analysis. The following logic is applied in determining whether to include the data in analysis:
  - If at least one of the eyes shows an abnormality in a particular parameter, the data for that specimen will be included in the statistical analysis even if the other eye is marked as "no data".
  - If the eyes are marked as "no data", or one eye is normal and the other eye is "no data" for a particular parameter the data for that specimen will not be included in the statistical analysis.

## Data QC

Image QC is typically performed during data collection to ensure high quality images are captured whilst eyes are dilated etc.

## Parameters and Metadata

**Left corneal thickness** TCPLA\_EYE\_066\_001 | v1.2

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

**VIP of right fundus** TCPLA\_EYE\_074\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

**Slit Lamp Equipment Model** TCPLA\_EYE\_032\_001 | v1.2

procedureMetadata

Req. Analysis: true      Req. Upload: false      Is Annotated: false

Options: SL130, BQ 900 LED/IM-900, SL 990, Micron III slit lamp extension, 30 SL-M, SL-15, SL-7E, SL30, S350, SL 139,

**Topical Anesthetic** TCPLA\_EYE\_044\_001 | v1.1

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Options:** Oxybuprocain, Phenylephrine hydrochloride, Mydriacyl, Atropine sulphate, No anesthesia, Hydrochloride, Atropine,

---

## Corneal Sclerization TCPLA\_EYE\_080\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** no data right eye, present both eyes, no data left eye, no data left eye, present right eye, present left eye, absent, present right eye, no data for both eyes, no data right eye, present left eye,

---

## Dilation Method TCPLA\_EYE\_043\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Options:** None, Tropicamide+Phenylephrin, Atropine sulphate, Atropine, Cyclopentolate hydrochloride+Phenylephrine hydrochloride, Cyclopentolate hydrochloride, Tropicamide, Phenylephrine hydrochloride,

---

## Min right eye lens density TCPLA\_EYE\_057\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** %

---

## Corneal ulcer TCPLA\_EYE\_085\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** no data for both eyes, absent, present right eye, no data left eye, present right eye, present left eye, present both eyes, no data left eye, no data right eye, no data right eye, present left eye,

---

## Iris/Pupil TCPLA\_EYE\_010\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** no data right eye, left eye abnormal, no data right eye, left eye abnormal, both eyes abnormal, no data left eye, no data left eye, right eye abnormal, no data for both eyes, normal, right eye abnormal,

---

## Experimenter ID TCPLA\_EYE\_036\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Iris transillumination TCPLA\_EYE\_082\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** left eye abnormal, both eyes abnormal, normal, no data left eye, no data left eye, right eye abnormal, no data right eye, right eye abnormal, no data right eye, left eye abnormal, no data for both eyes,

---

## Ophthalmoscope Lens Model TCPLA\_EYE\_089\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## VIP of right eye TCPLA\_EYE\_078\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Synechia TCPLA\_EYE\_019\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** present both eyes, no data right eye, no data for both eyes,  
no data left eye, present right eye, no data left eye, absent, present left eye,  
no data right eye, present left eye, present right eye,

---

## Pupil Light Response TCPLA\_EYE\_014\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** no data right eye, left eye abnormal, no data left eye, left eye abnormal,  
right eye abnormal, no data left eye, right eye abnormal, normal, no data right eye,  
both eyes abnormal, no data for both eyes,

---

## Scheimpflug description TCPLA\_EYE\_053\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

---

## Optical Coherence Tomography Equipment ID TCPLA\_EYE\_037\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Ophthalmoscope Observation TCPLA\_EYE\_029\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Corneal deposits TCPLA\_EYE\_081\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** present left eye, present both eyes, present right eye, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, absent, no data right eye, present left eye,

---

## Scheimpflug Equipment Manufacturer TCPLA\_EYE\_041\_001 | v1.4

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

**Options:** Oculus GmbH,

---

# Right corneal thickness TCPLA\_EYE\_060\_001 | v1.2

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

-----

# Images Slit Lamp TCPLA\_EYE\_051\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

-----

# Slit Lamp Equipment Manufacturer TCPLA\_EYE\_031\_001 | v1.2

procedureMetadata

Req. Analysis: true      Req. Upload: false      Is Annotated: false

Options: MuLe, Kowa, Haag-Streit, CSO, Zeiss, Topcon, Phoenix Research Labs,

-----

# Scheimpflug Equipment Model TCPLA\_EYE\_042\_001 | v1.4

procedureMetadata

Req. Analysis: true      Req. Upload: false      Is Annotated: false

**Options:** Pentacam,

---

**Max left eye lens density** TCPLA\_EYE\_055\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** %

---

**Retinal Blood Vessels Structure** TCPLA\_EYE\_025\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** left eye abnormal, both eyes abnormal, no data left eye, right eye abnormal, no data left eye, right eye abnormal, no data right eye, left eye abnormal, normal, no data right eye, no data for both eyes,

---

**General Anesthetic** TCPLA\_EYE\_045\_001 | v1.1

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Options:** No anesthesia, Ketamine+Medetomidine, Avertin, Isoflurane, Euthatal, Ketamine+Xylazine,

---

**Sheimpflug Lens description** TCPLA\_EYE\_052\_001 | v1.1

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

---

**Right anterior chamber depth** TCPLA\_EYE\_061\_001 | v1.2

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

**Unit Measured:** um

---

**Bulging eye** TCPLA\_EYE\_002\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

**Options:** no data right eye, present left eye, no data left eye, absent, present both eyes, no data left eye, present right eye, no data right eye, present left eye, present right eye, no data for both eyes,

---

## Fusion between cornea and lens TCPLA\_EYE\_018\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, present left eye, no data for both eyes,  
no data left eye, present right eye, no data right eye, present left eye, absent,  
present both eyes, present right eye, no data left eye,

---

## Pupil Position TCPLA\_EYE\_011\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** left eye abnormal, no data left eye, right eye abnormal, normal, no data right eye,  
no data right eye, left eye abnormal, right eye abnormal, no data left eye,  
no data for both eyes, both eyes abnormal,

---

## B-scan of left cornea and lens TCPLA\_EYE\_077\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

# Slit Lamp Equipment ID TCPLA\_EYE\_030\_001 | v1.2

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Lens TCPLA\_EYE\_016\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** no data left eye, no data for both eyes, left eye abnormal,  
no data left eye, right eye abnormal, no data right eye, no data right eye, left eye abnormal,  
right eye abnormal, normal, both eyes abnormal,

---

## Ophthalmoscope Equipment Manufacturer TCPLA\_EYE\_034\_001 | v1.2

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

**Options:** Heine / Volk, Kowa, Phoenix, Phoenix Research Labs, Haag-Streit, Heine,  
Keeler LTD, Karl Storz / Nikon,

---

## Eyelid closure TCPLA\_EYE\_005\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** left eye closed, no data left eye, no data left eye, right eye closed,  
no data for both eyes, no data right eye, right eye closed, normal, both eyes closed,  
no data right eye, left eye closed,

---

## Date OCT equipment last calibrated TCPLA\_EYE\_049\_001 | v1.1

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

## Lens Opacity TCPLA\_EYE\_017\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present both eyes, no data right eye, present left eye,  
no data left eye, present right eye, present left eye, no data for both eyes, no data left eye,  
absent, present right eye, no data right eye,

---

## Optical Coherence Tomography Equipment Manufacturer

TCPLA\_EYE\_038\_001 | v1.2

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** false

**Is Annotated:** false

**Options:** Bioptigen, Heidelberg Engineering,

---

**Vitreous** TCPLA\_EYE\_083\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** no data left eye, right eye abnormal, no data right eye, no data for both eyes, both eyes abnormal, no data left eye, right eye abnormal, left eye abnormal, no data right eye, left eye abnormal, normal,

---

**Ophthalmoscope Equipment Model** TCPLA\_EYE\_035\_001 | v1.2

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** false

**Is Annotated:** false

**Options:** Sigma 150K, Genesis-DF, Xenon Nova 175W light source + HOPKINS optic 1218AA /Nikon D5100 + 85 mm f/1.8 lens, Genesis-D, Micron III, OMEGA 180 / Superfield NC, Genesis, SL4 4AA, Omega 500 Unplugged, Omega 180 / 60D,

---

**Scheimpflug Equipment ID** TCPLA\_EYE\_040\_001 | v1.1

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

---

## VIP of left fundus TCPLA\_EYE\_075\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Corneal mineralization TCPLA\_EYE\_084\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** present right eye, no data right eye, absent, no data left eye, no data for both eyes, present left eye, no data right eye, present left eye, no data left eye, present right eye, present both eyes,

---

## Date Ophthalmoscope equipment last calibrated TCPLA\_EYE\_047\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Left total retinal thickness TCPLA\_EYE\_068\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Retinal Blood Vessels Pattern TCPLA\_EYE\_026\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data left eye, no data left eye, right eye abnormal, left eye abnormal, both eyes abnormal, right eye abnormal, no data for both eyes, normal, no data right eye, no data right eye, left eye abnormal,

---

## B-scan of right retina TCPLA\_EYE\_072\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Mean right eye lens density TCPLA\_EYE\_059\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

---

**Slit Lamp observation** TCPLA\_EYE\_028\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

**Persistence of hyaloid vascular system** TCPLA\_EYE\_027\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, present left eye, no data left eye, present right eye,  
no data for both eyes, no data left eye, present both eyes, absent, present left eye,  
no data right eye, present right eye,

---

**Images Ophthalmoscopy** TCPLA\_EYE\_050\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Eye TCPLA\_EYE\_001\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent both eyes, absent right eye, absent left eye,

---

## Optic Disc TCPLA\_EYE\_023\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: right eye abnormal, no data left eye, no data right eye, both eyes abnormal, normal, no data right eye, left eye abnormal, left eye abnormal, no data for both eyes, no data left eye, right eye abnormal,

---

## B-scan of left retina TCPLA\_EYE\_073\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Right total retinal thickness TCPLA\_EYE\_062\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

**Retina (combined)** TCPLA\_EYE\_092\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

---

**Retinal Blood Vessels** TCPLA\_EYE\_024\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

**Options:** left eye abnormal, no data right eye, both eyes abnormal,  
no data left eye, right eye abnormal, right eye abnormal, no data left eye,  
no data for both eyes, normal, no data right eye, left eye abnormal,

---

**VIP of left eye** TCPLA\_EYE\_079\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

**Eyelid morphology** TCPLA\_EYE\_004\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** normal, no data right eye, left eye abnormal, no data right eye, right eye abnormal, no data left eye, both eyes abnormal, left eye abnormal, no data for both eyes, no data left eye, right eye abnormal,

---

## Cornea TCPLA\_EYE\_007\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** left eye abnormal, no data for both eyes, no data left eye, right eye abnormal, no data right eye, both eyes abnormal, no data right eye, left eye abnormal, right eye abnormal, no data left eye, normal,

---

## Left eye diameter TCPLA\_EYE\_091\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Unit Measured:** mm

---

## Right posterior chamber depth TCPLA\_EYE\_065\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

## Min left eye lens density TCPLA\_EYE\_054\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

---

## Iris Pigmentation TCPLA\_EYE\_015\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** normal, no data for both eyes, both eyes abnormal, no data right eye, no data right eye, left eye abnormal, no data left eye, right eye abnormal, left eye abnormal, no data left eye, right eye abnormal,

---

## Pupil Dilation TCPLA\_EYE\_013\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data right eye, left eye dilated, no data right eye, no data for both eyes, no data left eye, right eye dilated, no data left eye, right eye dilated, normal, left eye dilated, both eyes dilated,

---

**Max right eye lens density** TCPLA\_EYE\_058\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

---

**Right inner nuclear layer** TCPLA\_EYE\_063\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

**Date Slit Lamp equipment last calibrated** TCPLA\_EYE\_046\_001 |

v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Corneal opacity TCPLA\_EYE\_008\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** present right eye, absent, present both eyes, no data left eye,  
no data for both eyes, no data left eye, present right eye, no data right eye,  
no data right eye, present left eye, present left eye,

---

## Optical Coherence Tomography Equipment Model TCPLA\_EY

E\_039\_001 | v1.2

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

**Options:** Spectralis, Envisu R2200, EnvisuTM R-Series SDOIS,

---

## Left anterior chamber depth TCPLA\_EYE\_067\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Unit Measured:** um

---

**Ophthalmoscope Equipment ID** TCPLA\_EYE\_033\_001 | v1.2

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

**Left vitreous humour thickness** TCPLA\_EYE\_088\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

---

**B-scan of right cornea and lens** TCPLA\_EYE\_076\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

**Narrow eye opening** TCPLA\_EYE\_006\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

**Options:** no data for both eyes, no data right eye, no data left eye, right eye abnormal, no data left eye, left eye abnormal, normal, right eye abnormal, both eyes abnormal, no data right eye, left eye abnormal,

---

**Left inner nuclear layer** TCPLA\_EYE\_069\_001 | v1.2

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: um

---

**Mean left eye lens density** TCPLA\_EYE\_056\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Unit Measured: %

---

**Lacrimation** TCPLA\_EYE\_086\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

**Options:** no data left eye, present right eye, no data left eye, present both eyes,  
no data right eye, present left eye, absent, present left eye, no data right eye,  
no data for both eyes, present right eye,

---

**Corneal vascularization** TCPLA\_EYE\_009\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

**Options:** present right eye, no data right eye, present left eye, absent, present left eye,  
no data right eye, no data for both eyes, no data left eye, no data left eye, present right eye,  
present both eyes,

---

**Right outer nuclear layer** TCPLA\_EYE\_064\_001 | v1.2

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

**Unit Measured:** um

---

**Date Scheimpflug equipment last calibrated** TCPLA\_EYE\_048\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

**Right vitreous humor thickness** TCPLA\_EYE\_087\_001 | v1.0  
simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

**Left posterior chamber depth** TCPLA\_EYE\_071\_001 | v1.2  
simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

---

**Right eye diameter** TCPLA\_EYE\_090\_001 | v1.0  
simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mm

---

# Pupil Shape TCPLA\_EYE\_012\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** right eye abnormal, no data for both eyes, no data left eye, right eye abnormal, normal, no data right eye, left eye abnormal, no data left eye, both eyes abnormal, no data right eye, left eye abnormal,

-----

# Eye Hemorrhage or Blood Presence TCPLA\_EYE\_003\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** no data left eye, present right eye, present right eye, no data right eye, no data left eye, no data right eye, present left eye, present left eye, present both eyes, no data for both eyes, absent,

-----

# Left outer nuclear layer TCPLA\_EYE\_070\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

-----