

MicroCT E14.5-E15.5 IMPC_EMO_001

Purpose

To assess morphological defects in E14.5-E15.5 embryos from lethal strains with MicroCT.

Experimental Design

- Set up timed matings with heterozygous mice
- Day 0 is defined as the midpoint of the prior dark cycle following the identification of a copulation plug.
- Collect embryos at E14.5 or E15.5 with 3 minimum homozygous and 1 littermate control
- Ensure embryos are bled as much as possible.
- Reconstruct the scans according to manufacturer's instructions.

Procedure

1. **Set up timed mating with heterozygous animals. Dissect and collect ≥ 3 homozygote embryos.**
2. **Process the embryos using the centre's staining and embedding protocol . Centre-specific methods are recorded within a document**
3. **Scan and generate reconstructions of embryos.**

Notes

Embryos may be processed for Histopathology.

Parameters and Metadata

Embryo reconstruction IMPC_EMO_001_001 | v1.2

seriesMediaParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Increments: Minimum 1

Experimenter ID IMPC_EMO_002_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Equipment ID IMPC_EMO_003_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Equipment Manufacturer IMPC_EMO_004_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Equipment Model IMPC_EMO_005_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Equipment Version IMPC_EMO_006_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Image Pixel Size IMPC_EMO_007_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: μm

E14.5-E15.5 MicroCT Preparation protocol IMPC_EMO_008_001 |

v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Object to Source IMPC_EMO_009_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: mm

Camera to Source IMPC_EMO_010_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: mm

Date equipment last calibrated IMPC_EMO_011_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Date of Sacrifice IMPC_EMO_012_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Date Preparation Finished IMPC_EMO_013_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Date of Scan IMPC_EMO_014_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Source Voltage IMPC_EMO_015_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: kV

Source Current IMPC_EMO_016_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: uA

Screenshots IMPC_EMO_017_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1
