

ALTERNATIVE - Ophthalmoscope ALTM-G-P_013_001

Parameters and Metadata

ALTERNATIVE - Eye ALTM-G-P_013_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Eye

Options: normal, abnormal,

ALTERNATIVE - Description ALTM-G-P_013_001_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Description

ALTERNATIVE - Left Retina ALTM-G-P_013_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Retina

Options: normal, abnormal, no data,

ALTERNATIVE - Left Retina Pigmentation ALTM-G-

P_013_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

ALTERNATIVE - Left Retina Structure ALTM-G-P_013_001_005 | v1.

0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Retina_Structure

Options: normal, degenerate, not defined, no data,

ALTERNATIVE - Left Optic Disc ALTM-G-P_013_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Optic_Disc

Options: normal, abnormal, not defined, no data,

ALTERNATIVE - Right Retina ALTM-G-P_013_001_007 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Retina

Options: normal, abnormal, no data,

ALTERNATIVE - Right Retina Pigmentation ALTM-G-

P_013_001_008 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

ALTERNATIVE - Right Retina Structure ALTM-G-P_013_001_009 |

v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Retina_Structure

Options: normal, degenerate, not defined, no data,

ALTERNATIVE - Right Optic Disc ALTM-G-P_013_001_010 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Optic_Disc

Options: normal, abnormal, not defined, no data,

ALTERNATIVE - Left Eye Blood Vessels ALTM-G-P_013_001_011 |

v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Eye_Blood_Vessels

Options: normal, abnormal, no data,

ALTERNATIVE - Left Eye Blood Vessels Structure ALTM-G-

P_013_001_012 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

ALTERNATIVE - Left Eye Blood Vessels Pattern ALTM-G-

P_013_001_013 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Left_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

ALTERNATIVE - Right Eye Blood Vessels ALTM-G-

P_013_001_014 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Eye_Blood_Vessels

Options: normal, abnormal, no data,

ALTERNATIVE - Right Eye Blood Vessels Structure ALTM-G-

P_013_001_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

ALTERNATIVE - Right Eye Blood Vessels Pattern ALTM-G-

P_013_001_016 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Right_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

ALTERNATIVE - Digital Image/s ALTM-G-P_013_001_017 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Digital_Image/s

ALTERNATIVE - Body Weight ALTM-G-P_013_001_018 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Body_Weight

ALTERNATIVE - Equipment name ALTM-G-P_013_001_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_name

ALTERNATIVE - Equipment manufacturer ALTM-G-

P_013_001_802 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_manufacturer

ALTERNATIVE - Equipment model ALTM-G-P_013_001_803 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_model

ALTERNATIVE - Topical agent(s) ALTM-G-P_013_001_804 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Topical_agent_s_

Options: atropine, tropicamide, atropine + neosynephrine, tropicamide + neosynephrine, tropicamide + phenylephrine,
