

ALTERNATIVE - Organ Weight ALTIMPC_OW T_001

Parameters and Metadata

ALTERNATIVE - Body weight ALTIMPC_OWT_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - body_weight

ALTERNATIVE - Spleen weight ALTIMPC_OWT_002_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - spleen_weight

ALTERNATIVE - Left kidney weight ALTIMPC_OWT_003_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - left_kidney_weight

ALTERNATIVE - Right kidney weight ALTIMPC_OWT_004_001 | v1.

0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - right_kidney_weight

ALTERNATIVE - Epididymis weight ALTIMPC_OWT_005_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - epididymis_weight

ALTERNATIVE - Seminal vesicle weight ALTIMPC_OWT_007_001

| v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - seminal_vesicle_weight

ALTERNATIVE - Spleen weight to body weight ratio ALTIMP

C_OWT_008_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - spleen_weight_to_body_weight_ratio

Derivation: $\text{div}(\text{'IMPC_OWT_002_001'}, \text{'IMPC_OWT_001_001'})$

ALTERNATIVE - Left kidney weight to body weight ratio AL

TIMPC_OWT_009_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - left_kidney_weight_to_body_weight_ratio

Derivation: $\text{div}(\text{'IMPC_OWT_003_001'}, \text{'IMPC_OWT_001_001'})$

ALTERNATIVE - Right kidney weight to body weight ratio

ALTIMPC_OWT_010_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - right_kidney_weight_to_body_weight_ratio

Derivation: div('IMPC_OWT_004_001', 'IMPC_OWT_001_001')

ALTERNATIVE - Epididymis weight to body weight ratio AL

TIMPC_OWT_012_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - epididymis_weight_to_body_weight_ratio

Derivation: div('IMPC_OWT_005_001', 'IMPC_OWT_001_001')

ALTERNATIVE - Seminal vesicle weight to body weight ratio

ALTIMPC_OWT_013_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - seminal_vesicle_weight_to_body_weight_ratio

Derivation: div('IMPC_OWT_007_001', 'IMPC_OWT_001_001')

ALTERNATIVE - Experimenter ID ALTIMPC_OWT_014_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - experimenter_id

ALTERNATIVE - Date of sacrifice ALTIMPC_OWT_015_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - date_of_sacrifice

ALTERNATIVE - Method of sacrifice ALTIMPC_OWT_016_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - method_of_sacrifice

ALTERNATIVE - Equipment ID ALTIMPC_OWT_017_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - equipment_id

ALTERNATIVE - Date equipment last calibrated ALTIMPC_OW

T_018_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - date_equipment_last_calibrated

ALTERNATIVE - Equipment manufacturer ALTIMPC_OWT_019_0

01 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - equipment_manufacturer

ALTERNATIVE - Equipment model ALTIMPC_OWT_020_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - equipment_model

ALTERNATIVE - Left testis weight ALTIMPC_OWT_021_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - left_testis_weight

ALTERNATIVE - Right testis weight ALTIMPC_OWT_022_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - right_testis_weight

ALTERNATIVE - Whole liver weight ALTIMPC_OWT_023_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - whole_liver_weight

ALTERNATIVE - Right lung all lobes weight ALTIMPC_OWT_024_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - right_lung_all_lobes_weight

ALTERNATIVE - Left lung weight ALTIMPC_OWT_025_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - left_lung_weight

ALTERNATIVE - Left testis weight to body weight ratio ALTI

MPC_OWT_026_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - left_testis_weight_to_body_weight_ratio

Derivation: $\text{div}(\text{'IMPC_OWT_021_001'}, \text{'IMPC_OWT_001_001'})$

ALTERNATIVE - Right testis weight to body weight ratio A

LTIMPC_OWT_027_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - right_testis_weight_to_body_weight_ratio

Derivation: $\text{div}(\text{'IMPC_OWT_022_001'}, \text{'IMPC_OWT_001_001'})$

ALTERNATIVE - Whole liver weight to body weight ratio AL

TIMPC_OWT_028_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - whole_liver_weight_to_body_weight_ratio

Derivation: div('IMPC_OWT_023_001', 'IMPC_OWT_001_001')

ALTERNATIVE - Right lung all lobes weight to body weight ratio ALTIMPC_OWT_029_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - right_lung_all_lobes_weight_to_body_weight_ratio

Derivation: div('IMPC_OWT_024_001', 'IMPC_OWT_001_001')

ALTERNATIVE - Left lung weight to body weight ratio ALTIMPC_OWT_030_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - left_lung_weight_to_body_weight_ratio

Derivation: div('IMPC_OWT_025_001', 'IMPC_OWT_001_001')

ALTERNATIVE - Whole lung weight ALTIMPC_OWT_031_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - whole_lung_weight
