

# Body Weight HRWLLA\_BWT\_001

## Purpose

The body weight test measures the weight of the mouse in a time series, allowing monitoring of its evolution; also, it is required in many other procedures.

## Equipment

- Laboratory Balance

## Procedure

- Place a cage containing mice on an operation table
- Remove a mouse from its cage. Determine and record its weight.
- Return the mouse to the cage or continue with a procedure.

## Notes

The body weight of each mouse is to be measured weekly. Optional additional weighs may also be recorded.

## Parameters and Metadata

### Experimenter ID HRWLLA\_BWT\_005\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

### General comments about the mouse HRWLLA\_BWT\_002\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Equipment ID HRWLLA\_BWT\_003\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

---

## Equipment manufacturer HRWLLA\_BWT\_004\_001 | v1.2

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Options:** Acculab, Sartorius, Mettler Toledo, Ohaus, Kern & Sohn GmbH, Scientech, Denver Instrument, Radwag,

---

## Body weight curve HRWLLA\_BWT\_008\_001 | v1.1

seriesMediaParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Unit Measured:** g

**Derivation:** unimplemented("")

**Increments:** Minimum 1

---

## Body weight HRWLLA\_BWT\_001\_001 | v1.3

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

---

## Date equipment last calibrated HRWLLA\_BWT\_006\_001 | v1.2

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Equipment model HRWLLA\_BWT\_007\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

**Options:** 440-47, 440-33, EW600-2M, Navigator 34120, 572-35, Entris 2201-1S, 440-47N, SL-3100D, VIC-511, 770-14, TP-202, PG3001-S, AV213C, AV2101, Adventurer Pro, AV212C, CPA3202S, EMB600-2, VIC-123, EMB500-1, EMB 200-2, XS802S, EMB220-1, PCB2000-1, AM100, BL310, BP6100, AB104-S, SI-2002, 201-10,

---