

Body Weight HAS_BWT_002

Purpose

The body weight test measures the weight of the mouse in a time series, allowing monitoring of its evolution; also, it is required in many other procedures.

Equipment

- Laboratory Balance

Procedure

- Place a cage containing mice on an operation table
- Remove a mouse from its cage. Determine and record its weight.
- Return the mouse to the cage or continue with a procedure.

Notes

The body weight of each mouse is to be measured weekly, from the 4th until the 16th week. Optional additional weighs may also be recorded.

Parameters and Metadata

Body weight HAS_BWT_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g

Description: body_weight

General comments about the mouse HAS_BWT_002_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: general_comments_about_the_mouse

Equipment ID HAS_BWT_003_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: equipment_name

Equipment manufacturer HAS_BWT_004_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: equipment_manufacturer

Options: Ohaus, Denver Instrument, Scientech, Mettler Toledo,

Experimenter ID HAS_BWT_005_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: experimenter_id

Date equipment last calibrated HAS_BWT_006_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Equipment model HAS_BWT_007_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: AV2101, SI-2002, SL-3100D, TP-202, PG3001-S,

Date of procedure HAS_BWT_008_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Procedural comments HAS_BWT_009_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

