

Slit Lamp ESLIM_014_001

Purpose

For viewing the eye structures in detail

Parameters and Metadata

Eye ESLIM_014_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: Eye

Options: normal, abnormal, left eye absent, right eye absent, bilaterally absent,

Description ESLIM_014_001_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: Description

Eye size ESLIM_014_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Eye_size

Options: normal, left eye decreased size, right eye decreased size, bilaterally decreased size, left eye increased size, right eye increased size, bilaterally increased size, not defined, no data,

Eye (+) blood presence ESLIM_014_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Eye__blood_presence

Options: absent, left eye present, right eye present, bilaterally present, not defined, no data,

Eyelid closure ESLIM_014_001_005 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Eyelid_closure

Options: normal, left eye closed, right eye closed, bilaterally closed, left eye closure incomplete, right eye closure incomplete, closure incomplete bilaterally,

Cornea ESLIM_014_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Cornea

Options: normal, left eye abnormal, right eye abnormal, bilaterally abnormal,

Cornea opacity ESLIM_014_001_007 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Cornea_opacity

Options: normal, left eye opaque, right eye opaque, bilaterally opaque, not defined, no data,

Cornea vascularisation ESLIM_014_001_008 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Cornea_vascularisation

Options: absent, left eye present, right eye present, bilaterally present, not defined, no data,

Lens ESLIM_014_001_009 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Lens

Options: normal, left eye abnormal, right eye abnormal, bilaterally abnormal,

Lens opacity ESLIM_014_001_010 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Lens_opacity

Options: normal, left eye opaque, right eye opaque, bilaterally opaque, not defined, no data,

Iris/pupil ESLIM_014_001_011 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: IrisPupil

Options: normal, right eye abnormal, left eye abnormal, bilaterally abnormal,

Iris/pupil position ESLIM_014_001_012 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: IrisPupil_position

Options: normal, left eye abnormal, right eye abnormal, bilaterally abnormal, not defined, no data,

Iris/pupil shape ESLIM_014_001_013 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: IrisPupil_shape

Options: normal, left eye abnormal, right eye abnormal, bilaterally abnormal, left eye dilated, right eye dilated, bilaterally dilated, not defined, no data,

Iris/pupil light response ESLIM_014_001_014 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: IrisPupil_light_response

Options: normal, left eye abnormal, right eye abnormal, bilaterally abnormal, no data,

Iris pigment ESLIM_014_001_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Iris_pigment

Options: normal, left eye abnormal, right eye abnormal, bilaterally abnormal, not defined, no data,

Digital image(s) ESLIM_014_001_016 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Digital_Image(s)

Equipment name ESLIM_014_001_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_name

Equipment manufacturer ESLIM_014_001_802 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_manufacturer

Equipment model ESLIM_014_001_803 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_model
