

Clickbox HAS_CLK_001

Parameters and Metadata

Date of procedure HAS_CLK_001_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

PIL number HAS_CLK_002_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

General comments about the mouse HAS_CLK_003_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: general_comments_about_the_mouse

Equipment ID HAS_CLK_004_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: equipment_name

Options: CB1, CB2, CB3, CB4, CB5, CB6,

Preyer Reflex HAS_CLK_005_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: preyer_reflex

Options: Normal (2), Reduced (1), Absent (0), Severely reduced (0.5),

Procedural comments HAS_CLK_006_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Body weight HAS_CLK_007_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: body_weight
